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
Fig 1. A: bright field image of HEK 293A cells, B: bright field image of HEK 293A cells with ROIs, C: SPR image of HEK 293A binding to hydroxystilbamidine. Bottom: distribution of K_D associated with cellular heterogeneity. Mean $K_D = 2.6 \mu\text{M}$.



Fig 2. A: bright field image of HEK 293A cells, B: bright field image of HEK 293A cells with ROIs, C: SPR image of HEK 293A binding to pentamidine. Bottom: distribution of K_D associated with cellular heterogeneity. Mean $K_D = 5.0 \mu\text{M}$.

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
Fig (a) and (b) VirD oscillators spotted on a SPR sensor chip and imaged with SPRM. Although the individual VirD oscillators cannot be resolved, the image intensity of each spot oscillates with the applied AC electrical field (c), providing the average oscillation amplitude of the VirD oscillators. The red line in (c) is the image intensity from the background regions, which was subtracted out in data analysis. (d) and (e) show the kinetics of a peptide ligand (SRIF14 MW=1560) binding to two different types of GPCRs (GABBR1 and NTSR1), where the solid lines are global fitting of the data to the first order kinetics. The extracted kinetic constants are shown on the graphs. Applied voltage: 0.3 V amplitude at 5 Hz. Buffer: PBS.

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1. Tontarawongsa S, Visitsattapongse S, Pechprasarn S. Performance Analysis of Non-Interferometry Based Surface Plasmon Resonance Microscopes. *Sensors (Basel)*, 2021,21(15).
2. Hassani H, Kreysing E. Noninvasive measurement of the refractive index of cell organelles using surface plasmon resonance microscopy. *Opt Lett*, 2019,44(6):1359-1362.
3. Kreysing E, Hassani H, Hampe N, et al. Nanometer-Resolved Mapping of Cell-Substrate Distances of Contracting Cardiomyocytes Using Surface Plasmon Resonance Microscopy. *ACS Nano*, 2018,12(9):8934-8942.
4. W. Jing et al, Simultaneous Quantification of Protein Binding Kinetics in Whole Cells with Surface Plasmon Resonance Imaging and Edge Deformation Tracking. *Membranes* 2020, 10, 247
5. H Yu, X Shan, S Wang, N Tao, Achieving high spatial resolution surface plasmon resonance microscopy with image reconstruction, *Anal. Chem.*, 2017, 89 (5), pp 2704-2707.
6. Y Wang, X Shan, H Wang, S Wang, N Tao, Plasmonic imaging of surface electrochemical reactions of single gold nanowires, *J. Am. Chem. Soc.*, 2017, 139 (4), pp 1376-1379.

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